

TITRATION LAB MANUAL

Experiment_ID : 23

Standardization of Acid by Titration with Standard Alkali

EXPERIMENT OBJECTIVE:

To standardize hydrochloric acid using standard sodium hydroxide solution

MATERIALS REQUIRED

Apparatus:

- * 2 Burettes - 50 mL each (Acid and Alkaline)
- * 3 Conical Flasks - 250 mL (Erlenmeyer)
- * 2 Beakers - 500 mL
- * 2 Pipettes - 20 mL Volumetric and Dropper
- * 1 Round-bottom Flask - 250 mL
- * 2 Burette Stands with clamps
- * 2 Glass Rods
- * 1 Funnel
- * 1 White Tile
- * Beaker Tongs

Chemicals:

- * Standard NaOH Solution - 0.1 M
- * HCl Solution - Unknown concentration
- * Phenolphthalein Indicator - 1% solution
- * Distilled Water

PROCEDURE - STEP BY STEP

STEP 1: Fill the Burette with NaOH

Instructions: Fill the **burette** with NaOH solution by selecting the **beaker** containing NaOH and using the **funnel** if required

Expected Outcome: Burette is filled with NaOH solution without air bubbles and ready for titration

STEP 2: Add HCl to the Conical Flask

Instructions: Add measured HCl solution into the **conical_flask** using the **pipette** from the **beaker**

Expected Outcome: Correct volume of HCl is present in the conical flask for titration

STEP 3: Add Indicator

Instructions: Add 2–3 drops of indicator into the **conical_flask** using the **dropper**

Expected Outcome: Solution remains colorless and ready to show endpoint color change

STEP 4: Perform the Titration

Instructions: Open the **burette** tap to add NaOH while swirling the **conical_flask** placed on the **white_tile**, and stop at the first permanent pale pink color

Expected Outcome: Solution changes from colorless to permanent pale pink indicating the endpoint

STEP 5: Record Final Burette Reading

Instructions: Record the final reading from the **burette** and note it in **observation**

Expected Outcome: Final reading is recorded to calculate the volume of NaOH used for titration

SAFETY PRECAUTIONS

SAFETY PRECAUTIONS AND GUIDELINES

- * Always wear safety glasses and lab coat when handling chemicals
- * Handle burettes and glassware carefully to avoid breakage
- * Never pipette by mouth - always use a rubber bulb
- * Ensure burette tips do not touch other glassware to avoid contamination
- * Add indicator drop by drop to avoid overdosing
- * Perform titration in well-ventilated area
- * Do not mix or store chemicals without proper labeling
- * In case of chemical splash, rinse immediately with water
- * Dispose of chemical waste in appropriate containers only
- * Keep MSDS (Material Safety Data Sheet) available
- * Know location of eye wash station and emergency equipment
- * Never leave heated equipment unattended
- * Wash hands thoroughly after handling chemicals
- * Report any spills or accidents to your instructor immediately

CALCULATIONS

Formula: $N_1V_1 = N_2V_2$

N_1 = Normality of HCl (to find)

V_1 = Volume of HCl used = 20 mL

N_2 = Normality of NaOH = 0.1 N

V_2 = Average volume of NaOH from concordant trials

To convert Normality to Molarity: $M = N / \text{number of } H^+/OH^-$

EXPECTED RESULTS

- * Typical HCl concentration: 0.1 - 0.2 M
- * Typical volume of NaOH used: 18-22 mL
- * Color change should be abrupt at endpoint
- * Concordant values within 0.20 mL difference

CONCLUSION

The titration experiment successfully standardizes HCl using standard NaOH. The concentration of HCl is determined using the formula $N_1V_1=N_2V_2$. Accurate results require careful technique, precision measurements, and adherence to all safety precautions.

----- END OF MANUAL -----

This manual covers all essential aspects of acid-base titration including procedures, expected outcomes, safety, and calculations.